

Built on the ACUVUE® MOIST Family—
the #1 prescribed daily disposable brand around the world¹

Fit Success Tips


1-DAY ACUVUE® MOIST Brand MULTIFOCAL Contact Lenses:
Helping you to continue meeting your patients' eye care needs as they age

1-DAY ACUVUE®
MOIST
BRAND CONTACT LENSES
MULTIFOCAL

Johnson & Johnson Vision Care have advanced the technology of multifocal lenses with 1-DAY ACUVUE® MOIST MULTIFOCAL

EYE-INSPIRED™ Design:

a unique aspheric centre-near design for the presbyopic eye



A consistent and predictable fit helps you continue to address your patients' needs as they age







Finding the best lens for the presbyopic eye












94% of patients were **successfully fit** with 2 pairs of lenses or less when using the Pocket Fitting Tips.*

Patient preparation:

1. Always start with a new subjective refraction, push the maximum plus for distance, and determine the spherical equivalent. Apply vertex distance correction if greater than +/-4.00D. Ensure that astigmatism is less than 1.00D.
2. Determine the dominant eye.
 - +1.00D blur test recommended
3. Determine the ADD based upon the patient's needs.
 - Focus on patient reading needs (desktop, tablet, mobile, etc)
4. Refer to the fit selection table based on the ADD to determine initial trial lenses.
 - Allow for 10 minutes of real-world exposure (outside of the exam room) before assessing visual performance
5. If your patient still has a specific near or distance need, refer to the enhancement tables.

To make your first lens the right lens, follow the selection table below:

For the best balance of distance and near vision		
ADD	EYE	LENS SELECTION
+0.75D to +1.25D	Dominant Eye	 LOW
	Non-dominant Eye	 LOW
+1.50D to +1.75D	Dominant Eye	 MID
	Non-dominant Eye	 MID
+2.00D to +2.50D	Dominant Eye	 MID
	Non-dominant Eye	 HIGH

For patients requiring an enhancement for distance vision or near vision	
ENHANCED DISTANCE VISION	ENHANCED NEAR VISION
1-DAY ACUVUE® MOIST BRAND SPHERE	 LOW
 LOW	 LOW +0.25D
 LOW	 MID
 MID	 MID +0.25D
 MID	 MID
 MID +0.25D	 HIGH +0.25D



Allow for 10 minutes of real-world exposure (outside of the exam room) before assessing visual performance.

*Up to 4 lenses total.

+: Add +0.25D to the distance power in the non-dominant eye.



1-DAY ACUVUE® MOIST MULTIFOCAL

Product Information

LENS DETAILS

Material	etafilcon A
Design	Aspheric centre-near design with INTUISIGHT™
Diameter	14.3 mm
Base curve	8.4 mm
Wetting agent	Embedded PVP/LACREON® Technology
UV blocker*	Class 2
Water content	58%
Dk†	21.4
Dk/t value‡	25.5
Centre thickness	0.084 mm for -3.00D (varies with power)
Visibility tint	Yes
Inside-out mark	Yes
Powers	-9.00D to +6.00D, 0.25D steps
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D

PVP=polyvinylpyrrolidone.

*All ACUVUE® Brand Contact Lenses have Class 1 or Class 2 UV blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area.

† Dk units: $\times 10^{-11}$ (cm²/sec) (ml O₂/ml x mm Hg). Dk/t units: 10-9 (cm/sec) (ml O₂/ml x mm Hg). All Dk values: Fatt units at 35°C, determined via polarographic method (boundary & edge corrected).

1. JJVC Data on file 2014. Internal analysis based on independent third party data, Nov 2013 to Oct 2014. 2. JJVC data on file 2014. 3. Sheardown H et al. Invest Ophthalmol Vis Sci 2006;47 [E-Abstract 2388]. 4. JJVC data on file 2005. 5. JJVC data on file 2007. 6. JJVC data on file 2009. [Post-hoc analysis, October 2009; adults aged 25–34 years (N=71); with 1-DAY ACUVUE® MOIST, 1% of patients experience frequent eye irritation, 14% of wearers occasionally experience eye irritation, while 83% seldom or never experienced eye irritation].

ACUVUE®, 1-DAY ACUVUE® MOIST, EYE-INSPIRED™, INTUISIGHT™, and LACREON® are trademarks of Johnson & Johnson Medical Ltd.
© Johnson & Johnson Medical Ltd. 2015. Johnson & Johnson Vision Care is part of Johnson & Johnson Medical Ltd.

1-DAY ACUVUE®
MOIST
BRAND CONTACT LENSES
MULTIFOCAL