

# EDUCATIONAL MOMENTS®



VISION CARE INSTITUTE

How to manage patients with CLPU (Contact Lens Peripheral Ulcer) or CLAIK (Contact Lens Associated Infiltrative Keratitis)

**INFLAMMATORY EVENT**

## WHAT YOU NEED TO KNOW

### Slit Lamp Viewing:


1. With fluorescein and yellow barrier filter. Optical section to assess depth
2. Medium/high magnification (16 - 25x)
3. Direct focal illumination

### Grading:

- Grade 0: None
- Grade 1: Trace
- Grade 2: Mild
- Grade 3: Moderate
- Grade 4: Severe

Location: superior, inferior, nasal, temporal

Central or peripheral



### Incidence:

- DW 1 - 5% (SiHs > to hydrogels)
- EW 1 - 13% (hydrogel), 3.3 - 5.4% (SiHs)

### Aetiology:

- Inflammatory reaction with focal excavation epithelium, infiltration and necrosis anterior stroma (Bowman's layer intact)
- Inflammatory response to exotoxins from gram positive bacteria (especially Staphylococcus spp)
- Due to: bacterial contamination, hypoxia, closed eye, tight lens, poor hygiene, solution toxicity, denatured lens deposits, mechanical trauma, lid margin disease
- Risk factors: high ametropia (>5D), younger age (15-25 years), case contamination, environmental influences, solution induced corneal staining (SICS)

### Symptoms:

- Up to 50% asymptomatic and unaware of complication
- Lens intolerance, foreign body sensation, photophobia, lacrimation, episodes of acute red eye
- Symptoms rapidly reduce after lens removal

### Signs:

- Moderate, localised hyperaemia
- Sterile infiltrate usually peripheral/mid-peripheral, small (0.2 - 2.0mm), single focal, circular, stains and diffuse infiltration (in acute phase) (can be multiple, although rare)
- No lid oedema; Anterior chamber reaction if severe

# EDUCATIONAL MOMENTS®



How to manage patients with CLPU (Contact Lens Peripheral Ulcer) or CLAIK (Contact Lens Associated Infiltrative Keratitis)

## WHAT YOU NEED TO RECOMMEND TO YOUR PATIENTS

### Recommendations:



- Close monitoring for 24 hours to ensure differential diagnosis from Microbial Keratitis (if central, >1mm and pain, treat suspiciously)
- Self limiting on removal
- Cease lens wear until epithelium intact over lesion (up to 14 days)
- Ocular lubricants to prevent lid rubbing and dilute bacterial toxins
- Oral analgesics to reduce discomfort
- Severe cases, acute red eye or no resolution with lens removal — refer for prophylactic antibiotic or combo-drug with mild anti-inflammatory (especially if infiltrates on visual axis)
- Eliminate bacterial source — lid hygiene, review hygiene
- Change to DW
- If recurrent, refit with RGP, increase lens replacement frequency, change care system
- If EW, consider limiting to 6 nights at a time
- Case hygiene (including rubbing & tissue wiping)
- For SICS, ensure a rub & rinse step, alter combination of SiH and MPS, switch to non-preserved solution or change to DD lens

### Prognosis:

- Good if visual axis not involved — single, spherical scar often remains, but can fade with time
- Certain subjects prone to recurrent inflammation (10 - 25%) — cease EW

### Differential Diagnosis:

Microbial keratitis (MK), marginal keratitis (Marg K), corneal dystrophies, corneal nerves, herpes simplex, stromal scar



NOTE: CLPU is also known as sterile corneal ulcer, staining infiltrate, sterile keratitis, sterile infiltrate

## HOW TO FIND OUT MORE


- ❖ Click [here](#) for a general refresher on slit lamp techniques
- ❖ Click [here](#) to watch our educational video on slit lamp examination using optical section

# EDUCATIONAL MOMENTS®



How to manage patients with CLPU (Contact Lens Peripheral Ulcer) or CLAIK (Contact Lens Associated Infiltrative Keratitis)

## PATIENT CASE STUDY



Patient DR is a 38-year-old car mechanic who has worn hydrogel daily wear lenses, seven days a week for 10 years.

He attends for a check-up for the first time in two years and initially reports no complaints about his lenses. On closer questioning he reveals that his right eye has been sore and red in the past few days. His fingers are stained with nicotine and his hands show signs of poor hygiene, with dirty fingernails. He is uncertain of the name of care solution he uses and no lens case has been brought to the appointment.

### Quiz:

**1. What grade would you give to this patient's CLPU?**

- A. Grade 1
- B. Grade 2
- C. Grade 3
- D. Grade 4

**2. What sign or symptom of CLPU is especially pertinent to its differential diagnosis?**

- A. Hyperaemia
- B. Mid peripheral/peripheral location
- C. Discomfort
- D. Lacrimation

**3. Which of the following signs and symptoms might determine whether you refer a patient with CLPU?**

- A. Central location
- B. >1mm in diameter
- C. Pain
- D. All of these

**4. Which of the following management options would you be most likely to choose first for this patient?**

- A. Discontinue lens wear and monitor closely for 24 hours before deciding on management
- B. Discontinue lens wear for two weeks before seeing again then decide on management
- C. Refer urgently for ophthalmological investigation
- D. Refit with DD lenses and advise on hand hygiene

Correct answers:

- 1: B Grade 2
- 2: B. CLPU is characterised by sterile infiltrates in a mid peripheral/peripheral location.
- 3: D. All of these signs and symptoms should be treated suspiciously to ensure differential diagnosis from MK.
- 4: A. CLPU is often self-limiting on lens removal but close monitoring for 24 hours is needed to exclude MK.

# EDUCATIONAL MOMENTS®



## How to manage patients with CLPU (Contact Lens Peripheral Ulcer) or CLAIK (Contact Lens Associated Infiltrative Keratitis )

### FURTHER READING/REFERENCES

Clinical Management Guidelines. Contact Lens-Associated Infiltrative Keratitis. College of Optometrists.

[CLICK HERE TO ACCESS](#)

Chalmers RL, Hickson-Curran SB, Keay L *et al.* Rates of adverse events with hydrogel and silicone hydrogel daily disposable lenses in a large postmarket surveillance registry: the TEMPO Registry. *Invest Ophthalmol Vis Sci* 2015;56:1 654-63.

[CLICK HERE TO ACCESS](#)

Szczotka-Flynn L, Jiang Y *et al.* Corneal inflammatory events with daily silicone hydrogel lens wear. *Optom Vis Sci* 2014;91:1 3-12.

[CLICK HERE TO ACCESS](#)

Szczotka-Flynn L and Chalmers RL. Incidence and epidemiologic associations of corneal infiltrates with silicone hydrogel contact lenses. *Eye Contact Lens* 2013;39:1 49-52.

[CLICK HERE TO ACCESS](#)

Chalmers RL, Keay L, McNally *et al.* Multicentre case-control study of the role of lens materials and care products on the development of corneal infiltrates. *Optom Vis Sci* 2012;89:3 316-25.

[CLICK HERE TO ACCESS](#)

Chalmers RL, Wagner H, Mitchell GL *et al.* Age and other risk factors for corneal infiltrative and inflammatory events in young soft contact lens wearers from the Contact Lens Assessment in Youth (CLAY) study. *Invest Ophthalmol Vis Sci* 2011;52:9 6690-6.

[CLICK HERE TO ACCESS](#)

Ozkan J, Mandathara P, Krishna P *et al.* Risk factors for corneal inflammatory and mechanical events with extended wear silicone hydrogel contact lenses. *Optom Vis Sci* 2010;87:11 847-53.

[CLICK HERE TO ACCESS](#)

Chalmers RL, Keay L, Long B *et al.* Risk factors for contact lens complications in US clinical practices. *Optom Vis Sci* 2010;87:10 725-35.

[CLICK HERE TO ACCESS](#)

Diec J, Carnt N, Tilia D *et al.* Prompt diagnosis and treatment of microbial keratitis in a daily wear lens. *Optom Vis Sci* 2009;86:7 E904-7.

[CLICK HERE TO ACCESS](#)

Radford CF, Minassian D, Dart JK *et al.* Risk factors for non-ulcerative contact lens complications in an ophthalmic accident and emergency department: a case-control study. *Ophthalmology* 2009;116:3 385-92.

[CLICK HERE TO ACCESS](#)

Sweeney DF and Naduvilath T. Are inflammatory events a marker for increased risk of microbial keratitis? *Eye Contact Lens* 2007;33:6 Pt 2 383-7.

[CLICK HERE TO ACCESS](#)

Stapleton F, Keay L, Jalbert I *et al.* The epidemiology of contact lens related infiltrates. *Optom Vis Sci* 2007;84:4 257-272.

[CLICK HERE TO ACCESS](#)

Efron N and Morgan PB. Can subtypes of contact lens associated corneal infiltrative events be clinically differentiated? *Cornea* 2006;25:5 540-4.

[CLICK HERE TO ACCESS](#)

Efron N and Morgan PB. Rethinking contact lens associated keratitis. *Clin Exp Optom* 2006;89:5 280-98.

[CLICK HERE TO ACCESS](#)

Morris J. Therapeutics in practice. The contact lens-related red eye: Part 2. *Optometry Today* 2004;10 September 32-37.

[CLICK HERE TO ACCESS](#)

Aasuri MK, Venkata M and Kumar VM. Differential diagnosis of microbial keratitis and contact lens-induced peripheral ulcer. *Eye Contact Lens* 2003;29:1 S60-2.

[CLICK HERE TO ACCESS](#)