

# NEW!

## EXPANDED PARAMETERS for a simpler fit process.



KEY MINUS POWERS				
Plano to -6.00 in 0.25D steps				
	-0.75	-1.25	-1.75	-2.25

PLUS POWERS AND HIGH MINUS POWERS			
-6.50 to -9.00 in 0.50D steps +0.25 to +4.00 in 0.25D steps			
	-0.75	-1.25	-1.75

OBLIQUES	30°	●	●	●			
	40°	●	●	●			
	50°	●	●	●			
	60°	●	●	●			
AGAINST THE RULE	70°	●	●	●	●	●	●
	80°	●	●	●	●	●	●
	90°	●	●	●	●	●	●
	100°	●	●	●	●	●	●
	110°	●	●	●	●	●	●
OBLIQUES	120°	●	●	●			
	130°	●	●	●			
	140°	●	●	●			
	150°	●	●	●			
WITH THE RULE	160°	●	●	●	●	●	●
	170°	●	●	●	●	●	●
	180°	●	●	●	●	●	●
	10°	●	●	●	●	●	●
	20°	●	●	●	●	●	●

● Current parameters    ● Expanded parameters

**EXPANDED PARAMETERS** mean  
**simplified experience** for your office<sup>1</sup>

ACUVUE® Brand 1-Day Contact Lenses for ASTIGMATISM have the **largest range of parameters** for astigmatism among daily disposable contact lenses<sup>1</sup>

1-DAY ACUVUE® MOIST Brand Contact Lenses for ASTIGMATISM now with **48% more parameters**<sup>1</sup>

Expanded parameters mean a **simplified experience** that helps you cover even more of your patients among daily disposable toric lenses<sup>1</sup>



Now covering

**87%**

**of astigmatic eyes**, expanding range from **1528** to **2260** parameters<sup>1</sup>

*1-DAY* ACUVUE®  
**MOIST**  
BRAND CONTACT LENSES  
FOR ASTIGMATISM

1. JJ&JV Data on file March 2017. Proportion astigmatism accommodated with ACUVUE® Brand CLs for ASTIGMATISM

ACUVUE® Brand Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available by visiting [jvisioncare.ae](http://jvisioncare.ae).

ACUVUE®, ACUVUE® MOIST are trademarks of Johnson & Johnson Middle East Inc. © Johnson & Johnson Middle East Inc. 2017.

171002024604916